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# INSTRUCTIONS FOR USE

## Anti-N (MNS2), monoclonal (mouse) OrthoClone

REF 690401A

### INTENDED USE

Monoclonal agglutinating Anti-N reagent is produced from cell culture supernatants of hybridoma cell lines. The cells are secreting an antibody, that reacts specifically with the corresponding antigen. The reagent is used to determine whether red blood cells possess or lack the corresponding blood group antigen. The reagent is intended to be used by qualified technical personnel only.

### PRINCIPLE OF PROCEDURE

The procedures used with this reagent are based on the principle of agglutination. Normal human erythrocytes, possessing the corresponding antigen, will agglutinate in the presence of the specific antibody directed towards the antigen.

### REAGENT

The listed reagent is produced by following clones:

Anti-N monoclonal, Mouse (clone 20H12; MN879)

This reagent contains <0.1% (w/v) sodium azide as preservative. Additionally the reagent is comprised of active antibody, sodium chloride, macromolecules and bovine albumin.

**CAUTION:** Please handle reagent with proper care: The reagent was prepared from supernatants of cell cultures. Handle this biological product as if capable of transmitting infectious agents.

The reagent contains sodium azide, which may be toxic and may react with lead or copper to form highly explosive salts. On disposal, flush with large quantities of water.

### STORAGE REQUIREMENTS

Store at +2 to +8 °C (unopened / opened), or at room temperature (+15 to +30 °C) while in use. Do not use reagent beyond its labeled expiration date.

### REMARKS

1. Strength of positive reactions also depends on the age of blood used.
2. It is recommended that each lot of reagent be tested with appropriate positive and negative controls.
3. Inappropriate storage impairs efficacy of the reagent.
4. Overcentrifugation or undercentrifugation may lead to false results.
5. Blood samples to be tested should be used as soon as possible. If a delay in testing occurs, samples should be stored at +2 to +8 °C. Blood drawn into sodium citrate or EDTA should be tested within 14 days.
6. The procedures identified below are for manual testing only. When using automated or semi-automated instruments, follow the procedures that are contained in the operator's manual provided by the device manufacturer. Laboratories must follow approved validation procedures to demonstrate compatibility of this product on automated systems.
7. For usage of this reagent all effective national laws, directives and guidelines must be followed..  
[In Germany especially: „Richtlinien zur Gewinnung von Blut und Blutbestandteilen und zur Anwendung von Blutprodukten (Hämotherapie)“.]

### REAGENT PREPARATION

There is no specific preparation of the reagent required. Use reagent directly from vials.

### PROCEDURE

Materials required but not provided:

at Slide Method

1. Glass slide
2. Pasteur pipette
3. Mixing stick

at Tube Centrifugation Method

1. Test tubes, 10 x 75 mm or 12 x 75 mm
2. Pipettes designed to deliver approximately 50 µL/100 µL
3. Centrifuge
4. Isotonic saline (with 0.85 - 0.9% sodium chloride)
5. Disposable pipette tips

### Test procedure

#### Slide Method

1. Use erythrocyte sediment or whole blood only.
2. Place one drop (approximately 50 µL) of reagent on a glass slide.
3. Using a Pasteur pipette add one drop erythrocyte sediment or whole blood (approximately 50 µL) to the glass slide.
4. Using a stick, mix the erythrocytes with reagent well and spread to a circle (diameter 2 cm).
5. By slightly rotating the slide, check for agglutination within 1 minute (reaction starts within seconds). Document the results. Unspecific reactions might appear due to drying of the reaction-formation or if the slide is heated.

#### Tube Centrifugation Method

1. Use a 2% to 5% suspension of red blood cells in isotonic saline only (cells washed one time or up to three times with isotonic saline).
2. Add 100 µL (alternative: one drop, approximately 50 µL) of reagent to each tube.
3. Add 100 µL (alternative: one drop, approximately 50 µL) of appropriate cell suspension to each tube.
4. Mix well by slightly shaking.
5. Incubate tube at room temperature (+15 to +30 °C) for 15 min.
6. Centrifugation of tube for 1 minute at 1.000 rpm (approximately 180 - 270 g).
7. Gently resuspend the red cells and within 3 minutes check macroscopically for agglutination. Document the results.

### INTERPRETATION OF RESULTS

"Slightly rotating/shaking " at Slide Method and at Tube Centrifugation Method.

**Positive result (+):** visible agglutination of erythrocytes indicates the presence of the corresponding antigen.

**Negative result (-):** no visible agglutination of erythrocytes indicates the absence of the corresponding antigen.

### LIMITATIONS OF THE PROCEDURE

1. The "Procedure" and "Interpretation of Results" sections must be followed closely to assure the accuracy of the test results.
2. No valid conclusion concerning the test result can be reached, if controls with uncertain or false results occur.
3. Enzyme treated erythrocytes may react unspecifically.
4. Addition of solutions containing bovine albumin or protein may lead to unspecific reactions.
5. Due to variability of antigen expression, reactivity of this reagent against certain phenotypes may produce a weaker reaction compared to control cells.
6. Red blood cells coated with alloantibodies or autoantibodies of the same or similar specificity as the reagent (i.e., cells that are positive in the direct antiglobulin-test [DAT]) may give weak reactions. In extreme cases, false-negative results may occur.
7. Anti-N may cause weak positive reactions with N-negative cells by the very similar biochemical basis of M and N antigen. This phenomenon occurs especially with commercial test cells and may be reduced by washing the cells.

### LEGEND OF SYMBOLS

	Store from - to	<b>REF</b>	Product Code		Expiration Date		Manufacturer as to 98/79/EU
<b>LOT</b>	Lot	<b>CLON</b>	Clone(s)	<b>IVD</b>	In vitro diagnostic medical device		EU CE-symbol

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## Anti-N (MNS2), monoclonal (mouse) OrthoClone

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