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# INSTRUCTIONS FOR USE

## Anti-H (A2), monoclonal IgM (mouse) OrthoClone

**REF** 690091A

### INTENDED USE

Monoclonal agglutinating Anti-H reagent is produced from cell culture supernatants of a hybridoma cell line. The cells secrete an antibody of IgM-type that reacts specifically with H-substance. The test serum is used to differentiate red blood cells into type A1 or A2. There is no free H-substance on A1 A1-bloods. A1B-, A1O- and especially BO-types have residual amounts of H-substance and show weak reactions with Anti-H. Strong reactions occur with A2-types or weaker A-types. As defined strongest reaction are seen with O-bloods. The reagent is intended to be used by qualified and technical personnel only.

### PRINCIPLE OF PROCEDURE

The procedure used with this reagent is based on the principle of agglutination. Normal human erythrocytes, possessing the corresponding antigen, will agglutinate in the presence of the specific antibody directed towards the antigen.

### REAGENT

The listed reagent contains antibodies of the following cell clone:

Anti- H monoclonal, mouse IgM (Clone: 10934C11).

The reagent contains <0.1% (w/v) sodium azide as preservative. Additionally the reagent is comprised of active antibody, sodium chloride, macromolecules and bovine albumin.

**CAUTION:** Please handle reagent with proper care: This reagent is prepared from supernatants of cell cultures. Handle this biological product as if capable of transmitting infectious agents. The reagent contains sodium azide, which may be toxic and may react with lead or copper to form highly explosive salts. On disposal, flush with large quantities of water.

### STORAGE REQUIREMENTS

Store at +2 to +8 °C (unopened / opened), or at room temperature (+15 to +30 °C) while in use. Do not use reagent beyond its labelled expiration date.

### REMARKS

1. Strength of positive reactions also depends on the age of blood used.
2. With each testing positive and negative controls should be performed.
3. Inappropriate storage impairs efficacy of the reagent.
4. Overcentrifugation or undercentrifugation may lead to false results.
5. Blood samples should be tested as soon as possible. If a delay in testing occurs, samples should be stored at +2 to +8 °C. Blood drawn into sodium citrate or EDTA should be tested within 14 days.
6. The procedure identified below is for manual testing only. When using automated or semi-automated instruments, follow the procedures that are contained in the operator's manual provided by the device manufacturer. Laboratories must follow approved validation procedures to demonstrate compatibility of this product on automated systems.
7. For usage of this reagent all effective national laws, directives and guidelines must be followed.  
[In Germany especially: „Richtlinien zur Gewinnung von Blut und Blutbestandteilen und zur Anwendung von Blutprodukten (Hämotherapie)“.]

### REAGENT PREPARATION

There is no specific preparation of the reagent required. Use reagent directly from the vial.

### PROCEDURE

Materials required but not provided:

At Tube Centrifugation Method

1. Test tubes, 10 x 75 mm or 12 x 75 mm
2. Pipettes designed to deliver 50 µL/100 µL
3. Centrifuge
4. Isotonic saline (0.85 – 0.9% sodium chloride)
5. Disposable pipette tips

### Test procedure

#### Tube Centrifugation Method

1. Use a 2% to 5% suspension of red blood cells in isotonic saline only (cells washed one time or up to three times with isotonic saline).
2. Add 100 µL (alternative: one drop, approximately 50 µL) of reagent to each tube.
3. Add 100 µL (alternative: one drop, approximately 50 µL) of appropriate cell suspension to each tube.
4. Mix well by shaking slightly.
5. Incubate tube at room temperature (+15 to +30 °C) for 15 min.
6. Centrifugation of tube for 1 minute at 2.000 rpm (approximately 800 – 1.000 g).
7. Resuspend the red cells by gently agitation and within 3 minutes check macroscopically for agglutination. Document the results.

### INTERPRETATION OF RESULTS

"Slightly shaking" at Tube Centrifugation Method.

**Positive result (+):** visible agglutination of erythrocytes indicates the presence of the corresponding antigen.

**Negative result (-):** no visible agglutination of erythrocytes indicates the absence of the corresponding antigen.

### LIMITATIONS OF THE PROCEDURE

1. The "Procedure" and "Interpretation of Results" sections must be followed closely to assure the accuracy of the test results.
2. No valid conclusion concerning the test result can be obtained, if controls with unclear or false results should occur.
3. Enzyme treated erythrocytes may react unspecifically.
4. Due to variability of antigen expression, reactivity of this reagent against certain phenotypes may produce a weaker reaction compared to control cells.
5. Red blood cells coated with alloantibodies or autoantibodies of the same or similar specificity as the reagent (i.e., cells that are positive in the direct antiglobulin-test [DAT]) may give weak reactions. In extreme cases, false-negative results may occur.

### SYMBOL - LEGENDE

 Store from - to	<b>REF</b> Product Code	 Expiration Date	 Manufacturer as to 98/79/EU
<b>LOT</b> Lot	<b>CLON</b> Clone(s)	<b>IVD</b> In vitro diagnostic medical device	 EU CE-symbol

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Per test in Provetta

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### Anti-H (A2), monoclonal IgM (mouse) OrthoClone

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