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### Anti-K (KEL1), monoclonal IgM (human) OrthoClone

**REF** 690461A

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 Antitoxin GmbH, Industriestraße 88, 69245 Bammental

730-22-0812 Version 012 / 30. Mrz 2015

**In vitro diagnostic Medical Device**

For Tube and Microplate Test and for combination with Ortho BioVue® Reverse Diluent Cassettes

# INSTRUCTIONS FOR USE

## Anti-K (KEL1), monoclonal IgM (human) OrthoClone

 690461A

### INTENDED USE

Monoclonal agglutinating Anti-K (KEL1) test serum is produced from cell culture supernatants of a heterohybridoma cell line. The cells secrete an antibody of IgM-type. The antibody is human protein. The test serum is used to determine whether red blood cells possess or lack the corresponding blood group antigen K. The procedures described are for the listed manual methods and application with Ortho BioVue® Reverse Diluent cassettes on the ORTHO VISION™ Analyzer and ORTHO AutoVue® Innova. The test serum is intended to be used by qualified technical personnel only.

### PRINCIPLE OF PROCEDURE

The procedures used with this reagent is based on the principle of agglutination. Normal human erythrocytes, possessing the corresponding antigen, will agglutinate in the presence of the specific antibody directed towards the antigen.

### REAGENT

The listed reagent is available with following clone:

**Anti-K** monoclonal, human IgM (clone: AEK-4)

This reagent contains <0.1% (w/v) sodium azide as preservative. Additionally the reagent is comprised of active antibody, sodium chloride, macromolecules and bovine albumin.

**CAUTION:** Please handle reagent with proper care: This reagent was prepared from supernatants of cell cultures. Handle as if capable of transmitting infectious agents. The reagent contains sodium azide, which may be toxic and may react with lead or copper to form highly explosive salts. On disposal, flush with large quantities of water.

### STORAGE REQUIREMENTS

Store at +2 to +8 °C (unopened / opened), or shortly at room temperature (+15 to +30 °C) while in use. If the devices ORTHO VISION™ Analyzer or ORTHO Auto Vue® Innova are used, the reagents may remain a maximum of 5 days on the device (+15 bis +30 °C). Do not use reagents beyond the stated expiration date!

### REMARKS

- Strength of positive reactions also depends on age of used blood.
- It is recommended that each lot of reagent be tested with appropriate positive and negative controls.
- Inappropriate storage impairs efficacy of the reagent.
- Overcentrifugation or undercentrifugation may lead to false results.
- Blood samples to be tested should be used as soon as possible. If a delay in testing occurs, samples should be stored at +2 to +8 °C. Blood drawn into sodium citrate or EDTA should be tested within 14 days.
- The procedures described are for the listed manual methods and analyzers only. When using other automated or semi-automated instruments, follow the device manufacturer's procedures in the operator's manual. Laboratories must follow approved validation procedures.
- For usage of this test serum all effective national laws, directives and guidelines must be followed. [In Germany especially: „Richtlinien zur Gewinnung von Blut und Blutbestandteilen und zur Anwendung von Blutprodukten (Hämotherapie)“.]

### REAGENT PREPARATION

There is no specific preparation of the reagent required. Use reagent directly from the vial.

### PROCEDURE

Material required but not provided:

#### Tube Centrifugation Method

- Test tubes, 10 x 75 mm or 12 x 75 mm
- Pipettes designed to deliver 50 µL/100 µL
- Centrifuge
- Isotonic saline (0,85 - 0,9% sodium chloride)
- Disposable pipette tips

#### Microplate Method

- Microplate with U-bottom, (optional: pre-treated)
- Centrifuge, suitable for centrifugation of microplates
- Microplate-carrier for centrifuge (option)
- Microplate shaker
- Microplatetest mirror (option)
- Pipette for 50 µL
- Isotonic saline (0.85 - 0.9 % sodium chloride)
- Disposable pipette tips

#### Application of BioVue® Cassettes

- Ortho BioVue® System centrifuge
- Isotonic saline (0,85 - 0,9 % sodium chloride)
- Mikroliter pipette for 10 µL, 40 µL and 50 µL
- Disposable pipette tips
- Ortho BioVue® System Work Rack
- Ortho BioVue® Reverse Diluent Cassettes

### Test procedure

#### Tube Centrifugation Method

- Use a 2% to 5% suspension of red blood cells in isotonic saline only (cells washed one time or up to three times with isotonic saline).
- Add 100 µL (alternative: one drop, approximately 50 µL) of reagent to each tube.
- Add 100 µL (alternative: one drop, approximately 50 µL) of appropriate cell suspension to each tube.
- Mix well by slightly shaking.
- Incubate tube at room temperature (+15 to +30 °C) for 15 min.
- Centrifugation of tube for 1 minute at 2.000 rpm (approximately 800-1.000 g).
- Gently resuspend the red cells and within 3 minutes check macroscopically for agglutination. Document the result.

### Microplate Method

1. Prepare a 2% to 5% suspension of red blood cells in isotonic saline. Cells can be washed up to three times with isotonic saline.
2. Add one drop (approximately 50 µL) of appropriate cell suspension into each well.
3. Add one drop (approximately 50 µL) of reagent into each well.
4. Mix microplate on a microplate shaker.  
**Note:** Recommended:  
(1) To mix: 10 to 30 seconds at medium strength.  
(2) To resuspend: 10 to 30 seconds at medium strength respectively with duration and strength, allowing a full resuspension of sediment without destruction of positive reactions.
5. Centrifugation of microplate.  
Recommended duration of centrifugation: 30 seconds at approximately 400 x g, or as long as necessary to produce best positive reactions with antigen positive cells, with used centrifuge and microplate, and still have an easy resuspension of antigen negative cells.  
**Note:** The used centrifugal force should be the minimum force to produce a sediment and a clear supernatant. Because of multitude offered centrifuges a general speed or duration for best results cannot be recommended. Each laboratory has to test and validate its equipment for speed and duration to obtain best results.
6. Resuspend the erythrocytes buttons
7. Position microplate for reading.
8. Read microplate and document the result.
9. Tests with negative or doubtful results have to be incubated for 5 to 10 minutes at room temperature (+15 to +30 °C).
10. Repeat steps 5 to 8 after incubation at room temperature.  
**Note:** All tests should be resuspended and evaluated directly after centrifugation.

### Manual Test Method in Ortho BioVue® Reverse Diluent Cassettes

The procedure described below is for manual Ortho BioVue® Cassettes testing only. Please refer to the instructions for use of the Ortho BioVue® Reverse Diluent Cassettes, sections STORAGE REQUIREMENT, PRECAUTIONS, SPECIMEN COLLECTION and PREPARATION, and INTERPRETATION OF RESULTS, as well as -where applicable- LIMITATIONS OF THE PROCEDURE.

1. For preparation of the cassettes, please refer to the instructions for use for Ortho BioVue® Reverse Diluent Cassettes, chapter „PROCEDURE“ point 1 and 2..
2. Prepare a 3 % to 5% erythrocyte suspension in isotonic sodium chloride solution. Optionally, the erythrocytes may be washed up to three times in isotonic sodium chloride solution.
3. Add 10 µL of a 3 to 5% suspension of test erythrocytes to the appropriate reaction chamber of the cassette.
4. Add 40 µL reagent to the appropriate reaction chambers of the cassette. The pipette tips shall not contact the rim of the reaction chambers. If this should adversely happen, please replace the tip with a new one, before you will fill the next reaction chamber.
5. Centrifuge the cassette in the Ortho BioVue® System centrifuge.  
**Note:** Centrifugation must be performed within 30 minutes after addition of specimens to the reaction chambers.
6. After completion of test, read each column from the cassette's front and rear sides for agglutination and/or hemolysis.
7. Protocol the final test result from the cassette side with the more clearly expressed positive reading.

### Test in Ortho BioVue® Reverse Diluent Cassettes on ORTHO VISION™ Analyzer or ORTHO AutoVue® Innova

Please follow the instructions in the user manual of the respective analyzer.

### INTERPRETATION OF RESULTS

"Slightly shaking" at Tube Centrifugation Method and using a microplate shaker at Microplate Method

**Positive result (+):** visible agglutination of erythrocytes indicates the presence of the corresponding antigen.

**Negative result (-):** no visible agglutination of erythrocytes indicates the absence of the corresponding antigen.

**Cassettes:** As for the manual method, for interpretation of results, please refer to the instructions for use of Ortho BioVue® Reverse Diluent Cassettes, chapter „INTERPRETATION OF RESULTS“. As for analyzers, reading and interpretation is effected automatically by each device.

### LIMITATIONS OF THE PROCEDURE

1. The "Procedure" and "Interpretation of Results" sections must be followed closely to assure the accuracy of the test results.
2. No valid conclusion concerning the test result can be reached, if controls with uncertain or false results occur.
3. Enzyme treated erythrocytes may react unspecifically.
4. Due to variability of antigen expression, reactivity of these reagents against certain phenotypes may produce a weaker reaction compared to control cells.
5. Red blood cells coated with alloantibodies or autoantibodies of the same or similar specificity as the reagent (i.e., cells that are positive in the direct antiglobulin-test [DAT]) may give weak reactions. In extreme cases, false-negative results may occur.
6. Red blood cells coated with antibodies (cells that are positive in the direct antiglobulin-test (DAT)) may give false-positive results in cassettes method. These cells react positive without reagent too.
7. Pay attention to all statements to limitations in inserts of used Ortho BioVue® Reverse Diluent Cassettes and the user manual of respective analyzer.

### PERFORMANCE

In compliance with Common Technical Specification (CTS Commission Decision of 03. February 2009) a performance evaluation was conducted for the product. Different samples (donor, patient, panel blood) was used and compared with other products. The determined sensitivity and specificity was as follows.

	Number of specimens tested	Specimens pos with Anti-K	Specificity	Specimens neg with Anti-K	Sensitivity
Tube / Mikroplate	310	196	100%	114	100%
Manual BioVue®	262	72	100%	190	99,47%
AutoVue®	102	16	100%	86	100%

### LEGEND OF SYMBOLS

	Store from - to		Product Code		Expiration Date		Manufacturer as to 98/79/EU
	Lot		Clone(s)		In vitro diagnostic medical device		EU CE-symbol

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