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INSTRUCTIONS FOR USE

Anti-Kell (KEL1), Coombsreactive, polyclonal (human)

REF 690482A

INTENDED USE

Polyclonal Coombsreactive Anti-K (KEL1) reagent is produced from human plasma that contains a specific antibody of IgG-type, which reacts exclusively with the corresponding antigen. The reagent is used to determine whether red blood cells possess or lack the corresponding blood group antigen K. The reagent is intended to be used by qualified and technical personnel only.

PRINCIPLE OF PROCEDURE

The procedure used with this reagent is based on the principle of agglutination. Normal human erythrocytes, possessing the antigen, will be recognized and coated with the corresponding specific antibody and then the cells will be agglutinated by a secondary antibody, which reacts with human IgG-molecules.

REAGENT

The listed reagent is available in following formulation:

Anti-K (KEL1) Coombsreactive (polyclonal, human).

The reagent contains <0.1% (w/v) sodium azide as preservative. Beside the active antibody and the human serum, the reagent contains sodium chloride, macromolecules and bovine albumin.

CAUTION: Please handle reagent with proper care: This reagent is prepared from human plasma. The raw materials for this product are tested for HBsAg, HIV- and HCV antibodies and found to be negative. Handle as if capable of transmitting infectious agents. The reagent contains sodium azide, which may be toxic and may react with lead or copper to form highly explosive salts. On disposal, flush with large quantities of water.

STORAGE REQUIREMENTS

Store at +2 to +8 °C (unopened / opened), or at room temperature (+15 to +30 °C) while in use. Do not use reagent beyond its labelled expiration date!

REMARKS

1. Strength of positive reactions also depends on age of blood used.
2. It is recommended that each lot of reagent be tested with appropriate positive and negative controls.
3. Weak turbidity of the reagent does not affect its reactivity.
4. Inappropriate storage impairs efficacy of the reagent.
5. Overcentrifugation or undercentrifugation may lead to false results.
6. Blood samples to be tested should be used as soon as possible. If a delay in testing occurs, samples should be stored at +2 to +8 °C. Blood drawn into sodium citrate or EDTA should be tested within 14 days.
7. The procedure identified below is for manual testing only. When using automated or semi-automated instruments, follow the procedures that are contained in the operator's manual provided by the device manufacturer. Laboratories must follow approved validation procedures to demonstrate compatibility of this product on automated systems.
8. For usage of this reagent all effective national laws, directives and guidelines must be followed.
[In Germany especially: „Richtlinien zur Gewinnung von Blut und Blutbestandteilen und zur Anwendung von Blutprodukten (Hämotherapie)“.]

REAGENT PREPARATION

There is no specific preparation of the reagent required. Use reagent directly from the vials.

PROCEDURE

Material required but not provided:

At Tube Centrifugation Method

1. Test tubes, 10 x 75 mm or 12 x 75 mm
2. Pipettes designed to deliver 50 µL/100 µL
3. Centrifuge
4. Isotonic saline (0.85 – 0.9% sodium chloride)
5. Anti-Human-Globulin-serum (Coombs-serum)
6. Disposable pipette tips

Test procedure

Tube Centrifugation Method

1. Use a 2% to 5% suspension of red blood cells in isotonic saline only (cells washed one time or up to three times with isotonic saline).
2. Add 100 µL (alternative: one drop, approximately 50 µL) of reagent to each tube.
3. Add 100 µL (alternative: one drop, approximately 50 µL) of appropriate cell suspension to each tube.
4. Mix well by slightly shaking.
5. Incubate tubes at +37 °C for 30 min.
6. Wash red cells 3 times with (cold) isotonic saline.
7. Add 100 µL Anti-Human-Globulin-serum (Coombs-serum) to each tube.
8. Centrifugation of tube for 1 min at 1.000 rpm (approximately 180 - 270 g).
9. Gently resuspend the red cells and within 3 minutes check macroscopically for agglutination. Document the results.

INTERPRETATION OF RESULTS

"Slightly shaking" at Tube Centrifugation Method

Positive result (+): visible agglutination of erythrocytes indicates the presence of the corresponding antigen.

Negative result (-): no visible agglutination of erythrocytes indicates the absence of the corresponding antigen.

LIMITATIONS OF THE PROCEDURE

1. The "Procedure" and "Interpretation of Results" sections must be followed closely to assure the accuracy of the test results.
2. Enzyme treated erythrocytes may react unspecific.
3. Due to variability of antigen expression, reactivity of this reagent against certain phenotypes may produce a weaker reaction compared to control cells.
4. Red blood cells coated with alloantibodies or autoantibodies of the same or similar specificity as the reagent (i.e., cells that are positive in the direct antiglobulin-test [DAT]) are not suitable for this test procedure.
5. The described procedure for testing is for usage with Coombs-serum produced by Antitoxin GmbH. On principle Coombs-sera of other manufacturers can be used, but laboratories have to follow the procedures of this manufacturer. Also laboratories must follow approved validation procedures to demonstrate compatibility of this product.
6. No valid conclusion concerning the test result can be obtained, if controls with unclear or false results should occur.

PERFORMANCE

In compliance with Common Technical Specifications (CTS Commission Decision of 03. February 2009) a performance evaluation was conducted. Different samples (donor, patient, panel blood) was used and compared with other products.

Product	Positive	False negative	Sensitivity	Negative	False positive	Specificity
Anti-Kell Coombsreactive	73	0	100%	511	0	100%

SYMBOL - LEGENDE

 Store from - to	REF Product Code	 Expiration Date	 Hersteller nach 98/79/EU
LOT Lot	CLON Clone(s)	IVD In vitro diagnostic medical device	 EU CE-symbol

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Anti-Kell (KEL1), coombsreactive, polyclonal (human)

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