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# INSTRUCTIONS FOR USE

**Anti-Lu<sup>a</sup>, Anti-Lu<sup>b</sup>; coombsreactive, polyclonal (human)**  
**Anti-Xg<sup>a</sup>; coombsreactive, polyclonal (human)**  
**Anti-Kp<sup>a</sup>; coombsreactive, polyclonal (human)**

**REF** 690612A und 690622A**REF** 690632A**REF** 690652A

## INTENDED USE

Polyclonal Coombsreactive Anti-Lu<sup>a</sup>, Anti-Lu<sup>b</sup>, Anti-Xg<sup>a</sup> and Anti-Kp<sup>a</sup> reagents are produced from human plasma containing a specific antibody of IgG-type, which reacts exclusively with the corresponding antigen. The reagents are used to determine whether red blood cells possess or lack the corresponding blood group antigens Lu<sup>a</sup>, Lu<sup>b</sup>, Xg<sup>a</sup> and Kp<sup>a</sup>. The reagents are intended to be used by qualified technical personnel only.

## PRINCIPLE OF PROCEDURE

The procedure used with these reagents is based on the principle of agglutination. Normal human erythrocytes, possessing one of these antigens, will be recognized and coated with the corresponding specific antibody and then the cells will be agglutinated by a secondary antibody, that reacts with human IgG-molecules.

## REAGENT

The listed reagents are comprised of, with differences in specificity for the antibody only:

Anti-Lu<sup>a</sup>, Anti-Lu<sup>b</sup>, Anti-Xg<sup>a</sup>, Anti-Kp<sup>a</sup> Coombsreactive, polyclonal (human)

All reagents contain <0.1% (w/v) sodium azide as preservative. Besides the active antibody and human serum, the reagents contain sodium chloride, macromolecules and bovine albumin.

**CAUTION:** Please handle reagents with proper care. These reagents were prepared from human plasma. The raw materials for these products are tested for HBsAg, HIV- and HCV antibodies and found to be negative. Handle as if capable of transmitting infectious agents. The reagents contain sodium azide, which may be toxic and may react with lead or copper to form highly explosive salts. On disposal, flush with large quantities of water.

## STORAGE REQUIREMENTS

Store at +2 to +8 °C (unopened / opened), or at room temperature (+15 to +30 °C) while in use. Do not use reagents beyond their labelled expiration date!

## REMARKS

1. Strength of positive reactions also depends on age of blood used.
2. It is recommended that each lot of reagents be tested with appropriate positive and negative controls.
3. Weak turbidity of the reagents does not affect its reactivity.
4. Inappropriate storage impairs efficacy of the reagents.
5. Overcentrifugation or undercentrifugation may lead to false results.
6. Blood samples to be tested should be used as soon as possible. If a delay in testing occurs, samples should be stored at +2 to +8 °C. Blood drawn into sodium citrate or EDTA should be tested within 14 days.
7. The procedure identified below is for manual testing only. When using automated or semi-automated instruments, follow the procedures that are contained in the operator's manual provided by the device manufacturer. Laboratories must follow approved validation procedures to demonstrate compatibility of this product on automated systems.
8. For usage of these reagents all effective national laws, directives and guidelines must be followed.  
[In Germany especially: „Richtlinien zur Gewinnung von Blut und Blutbestandteilen und zur Anwendung von Blutprodukten (Hämotherapie)“.]

## REAGENT PREPARATION

There is no specific preparation of the reagents required. Use reagents directly from the vials.

## PROCEDURE

Material required but not provided at:

### Tube Centrifugation Method

1. Test tubes, 10 x 75 mm or 12 x 75 mm
2. Pipettes designed to deliver 50 µL/100 µL
3. Centrifuge
4. Isotonic saline (0.85 – 0.9% sodium chloride)
5. Anti-Human-Globulin-serum (Coombs-serum)
6. Disposable pipette tips

## Test procedure

### Tube Centrifugation Method

1. Use a 2% to 5% suspension of red blood cells in isotonic saline only (cells washed one time or up to three times with isotonic saline).
2. Add 100 µL (alternative: one drop, approximately 50 µL) of appropriate reagent to each tube.
3. Add 100 µL (alternative: one drop, approximately 50 µL) of appropriate cell suspension to each tube.
4. Mix well by slightly shaking.
5. Incubate tubes at +37 °C for 30 min.
6. Wash red cells 3 times with (cold) isotonic saline.
7. Add 100 µL Anti-Human-Globulin-serum (Coombs-serum) to each tube.
8. Centrifugation of tubes for 1 min at 1.000 rpm (approximately 180 - 270 g).
9. Gently resuspend the red cells and within 3 minutes check macroscopically for agglutination. Document the results.

## INTERPRETATION OF RESULTS

"Slightly shaking" at Tube Centrifugation Method

**Positive result (+):** visible agglutination of erythrocytes indicates the presence of the corresponding antigen.

**Negative result (-):** no visible agglutination of erythrocytes indicates the absence of the corresponding antigen.

## LIMITATIONS OF THE PROCEDURE

1. The "Procedure" and "Interpretation of Results" sections must be followed closely to assure the accuracy of the test results.
2. Enzyme treated erythrocytes may react unspecific.
3. Due to variability of antigen expression, reactivity of these reagents against certain phenotypes may produce a weaker reaction compared to control cells.
4. Red blood cells coated with alloantibodies or autoantibodies of the same or similar specificity as the reagent (i.e., cells that are positive in the direct antiglobulin-test [DAT]) are not suitable for this test procedure.
5. The described procedure for testing is for usage with Coombs-serum produced by Antitoxin GmbH. On principle Coombs-sera of other manufacturers can be used, but laboratories have to follow the procedures of this manufacturer. Also laboratories must follow approved validation procedures to demonstrate compatibility of this product.
6. No valid conclusion concerning the test result can be obtained, if controls with unclear or false results should occur.

## SYMBOL - LEGENDE

	Store from - to	<b>REF</b>	Product Code		Expiration Date		Hersteller nach 98/79/EU
<b>LOT</b>	Lot	<b>CLON</b>	Clone(s)	<b>IVD</b>	In vitro diagnostic medical device		EU CE-symbol

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**Anti-Lu<sup>a</sup>, Anti-Lu<sup>b</sup>; coombsreactive, polyclonal (human)**  
**Anti-Xg<sup>a</sup>; coombsreactive, polyclonal (human)**  
**Anti-Kp<sup>a</sup>; coombsreactive, polyclonal (human)**

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