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Anti-D (RH1), P3x61, monoclonal IgM (human)

REF 690121A

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REF 690121A

INTENDED USE

Monoclonal agglutinating Anti-D reagent is produced from cell culture supernatants of hetero-hybridoma cell lines. The cells secrete an antibody of IgM-type that reacts specifically with the corresponding antigen. The antibody is human protein. The reagent is used to determine whether human red blood cells possess or lack the corresponding blood group antigen D. This reagent reacts negatively with erythrocytes of category D^{VI}. The reagent is intended to be used by qualified and technical personnel only.

PRINCIPLE OF PROCEDURE

The procedures used with this reagent is based on the principle of agglutination. Normal human erythrocytes, possessing the corresponding antigen, will agglutinate in the presence of the specific antibody directed towards the antigen.

REAGENT

The listed reagent contains antibodies of the following cell clone:

Anti-D monoclonal, human IgM (clone: P3X61).

This reagent contains <0.1% (w/v) sodium azide as preservative. Additionally the reagent is comprised of active antibody, sodium chloride, macromolecules and bovine albumin.

CAUTION: Please handle reagent with proper care: This test serum was prepared from supernatants of cell cultures. As biological product it should be looked upon as potentially infectious because of never complete exclusion of danger through excitants of disease. The test serum contains sodium azide that may be toxic and may react with lead or copper to form highly explosive salts. On disposal, flush with large quantities of water. Found from the above reasons test serum should be handled with proper care.

STORAGE REQUIREMENTS

Store at +2 to +8 °C (unopened / opened), or at room temperature (+15 to +30 °C) while in use. Do not use reagent beyond its labelled expiration date!

REMARKS

- Strength of positive reactions also depends on age of used blood.
- It is recommended that each lot of reagent be tested with appropriate positive and negative controls.
- Inappropriate storage impairs efficacy of the reagent.
- Overcentrifugation or undercentrifugation may lead to false results.
- Blood samples to be tested should be used as soon as possible. If a delay in testing occurs, samples should be stored at +2 to +8 °C. Blood drawn into sodium citrate or EDTA should be tested within 14 days.
- The procedures identified below are for manual testing only. When using automated or semi-automated instruments, follow the procedures that are contained in the operator's manual provided by the device manufacturer. Laboratories must follow approved validation procedures to demonstrate compatibility of this product on automated systems.
- For usage of this reagent all effective national laws, directives and guidelines must be followed.
[In Germany especially: „Richtlinien zur Gewinnung von Blut und Blutbestandteilen und zur Anwendung von Blutprodukten (Hämotherapie)“.]
- This test serum reacts negatively with erythrocytes of category D^{VI}.

REAGENT PREPARATION

There is no specific preparation of the reagent required. Use reagent directly from the vials.

PROCEDURE

Material required but not provided:

Slide Method

- Glass slide
- Pasteur pipette
- Mixing stick

Tube Centrifugation Method

- Test tubes, 10 x 75 mm or 12 x 75 mm
- Pipettes designed to deliver 50 µL/100 µL
- Centrifuge
- Isotonic saline (0,85 - 0,9% Natriumchlorid)
- Disposable pipette tips

Test procedure

Slide Method

- Use erythrocyte sediment or whole blood only.
- Place one drop (approximately 50 µL) of reagent on a glass slide.
- Using a Pasteur pipette add one drop erythrocyte sediment or whole blood (approximately 50 µL) to the glass slide.
- Using a stick, mix the erythrocytes with reagent well and spread to a circle (diameter 2 cm).
- By slightly rotating the slide, check for agglutination within 1 minute (reaction starts within seconds). Document the results. Unspecific reactions might appear due to drying of the reaction-formation or if the slide is heated.

Tube Centrifugation Method

- Use a 2% to 5% suspension of red blood cells in isotonic saline only (cells washed one time or up to three times with isotonic saline).
- Add 100 µL (alternative: one drop, approximately 50 µL) of reagent to each tube.
- Add 100 µL (alternative: one drop, approximately 50 µL) of appropriate cell suspension to each tube.
- Mix well by slightly shaking.
- Incubate tube at room temperature (+15 to +30 °C) for 1-15 min.
- Centrifugation of tube for 1 minute at 2.000 rpm (approximately 800 - 1.000 g).
- Gently resuspend the red cells and within 3 minutes check macroscopically for agglutination. Document the result.

INTERPRETATION OF RESULTS

"Slightly rotating / shaking" at Slide Method and at Tube Centrifugation Method.

Positive result (+): visible agglutination of erythrocytes indicates the presence of the corresponding antigen.

Negative result (-): no visible agglutination of erythrocytes indicates the absence of the corresponding antigen.

LIMITATIONS OF THE PROCEDURE

- The "Procedure" and "Interpretation of Results" sections must be followed closely to assure the accuracy of the test results.
- No valid conclusion concerning the test result can be reached, if controls with uncertain or false results occur.
- Enzyme treated erythrocytes may react unspecifically.
- Due to variability of antigen expression, reactivity of this reagent against certain phenotypes may produce a weaker reaction compared to control cells.
- Red blood cells coated with alloantibodies or autoantibodies of the same or similar specificity as the reagent (i.e., cells that are positive in the direct antiglobulin-test [DAT]) may give weak reactions. In extreme cases, false-negative results may occur.
- With the slide method will most of weak D-antigens (D^{weak}) and categories not be recognized.
- This test serum reacts negatively with erythrocytes of category D^{VI}.

PERFORMANCE

In compliance with Common Technical Specifications (CTS Commission Decision of 03. February 2009) a performance evaluation was conducted. Different samples (donor, patient, panel blood) was used and compared with other products. The calculated diagnostic values found for these studies:

Product	Positive	False negative	Sensitivity	Negative	False positive	Specificity
Anti-D (clone: P3x61)	1986	0*	100%*	517	0	100%

* for Anti-D (clone: P3x61) in slide method 8 D^{weak} samples were tested false negative, so in this single method sensitivity was 99,56% only.

LEGEND OF SYMBOLS

 Store from - to	REF	Product Code	 Expiration Date	 Manufacturer as to 98/79/EU
LOT Lot	CLON	Clone(s)	IVD In vitro diagnostic medical device	 EU CE-symbol

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