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## Anti-CDE (RH2,1,3), monoclonal IgM + IgG (human) OrthoClone

**REF** 690141A

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 Antitoxin GmbH, Industriestraße 88, 69245 Bammental

730-22-6306 Version 006 / 01. Aug. 2014

**In vitro diagnostic Medical Device**  
For Slide and Tube Test

## INSTRUCTIONS FOR USE

### Anti-CDE (RH2,1,3), monoclonal IgM + IgG (human) OrthoClone

 690141A

#### INTENDED USE

Monoclonal agglutinating Anti-CDE test serum is produced from cell culture supernatants of heterohybridoma cellines. The cells are secreting an antibody of IgM- or IgG-type. Each antibody is human protein. The test serum is used to determine whether red blood cells possess or lack the blood group antigens C, D and E. The test serum is intended to be used by qualified and technical personnel only.

#### PRINCIPLE OF PROCEDURE

The procedures used with this test serum are based on the principle of agglutination. Normal human erythrocytes, possessing one of the corresponding antigens, will agglutinate in the presence of the specific antibody directed towards the antigen.

#### REAGENT

The listed test serum is produced by following cell clones:

- Anti-C (clone: MS-24 (IgM))
- Anti-D (clones: MS-201 (IgM), MS-26 (IgG))
- Anti-E (clone: MS-80 (IgM))

The test serum contains <0.1% (w/v) sodium azide as preservative. Additionally the test serum is comprised of active antibody, sodium chloride, macromolecules and bovine albumin.

**CAUTION:** This test serum was prepared from supernatants of cell cultures. As biological product it should be looked upon as potentially infectious because of never complete exclusion of danger through excitants of disease. The test serum contains sodium azide that may be toxic and may react with lead or copper to form highly explosive salts. On disposal, flush with large quantities of water. Found from the above reasons test serum should be handled with proper care.

#### STORAGE REQUIREMENTS

Store at +2 to +8°C (unopened/opened), or at room temperature (+15 to +30°C) while in use. Do not use reagent beyond their labeled expiry date.

#### REMARKS

- Strength of positive reactions also depends on age of used blood
- It is recommended that each lot of test serum should be tested with appropriate positive and negative controls.
- Inappropriate storage impairs efficacy of the test serum.
- Centrifugation highly different from appointed relative centrifugal force may lead to false results
- Blood samples to be tested should be used as soon as possible. If a delay in testing occurs, samples should be stored at +2 to +8°C. Blood drawn into sodium citrate or EDTA should be tested within 14 days.
- The procedures identified below are for manual testing only. When using automated or semi-automated instruments, follow the procedures that are contained in the operator's manual provided by the device manufacturer. Laboratories must follow approved validation procedures to demonstrate compatibility of this product on automated systems.
- For usage of this test serum all effective national laws, directives and guidelines must be followed.  
[in Germany especially the „Richtlinien zur Gewinnung von Blut und Blutbestandteilen und zur Anwendung von Blutprodukten (Hämotherapie)“]

#### REAGENT PREPARATION

There is no specific preparation of the test serum required. Use test serum directly from the vial.

#### PROCEDURE

Material required but not provided:

##### Slide Method

- Glass slide
- Pasteur pipette
- Mixing stick

##### Tube Centrifugation Method

- Test tubes, 10 x 75 mm or 12 x 75 mm
- Pipettes designed to deliver 100 µL, 50 µl
- Centrifuge
- Isotonic saline (0.85 – 0.9% sodium chloride)
- Disposable pipette tips
- Anti-Human-Globulin-Serum (Coombs-serum)

#### Test procedure

##### Slide Method

- Use erythrocyte sediment only.
- Place one drop (approximately 50 µL) of test serum on a glass slide.
- Using a Pasteur pipette add one drop erythrocyte sediment (approximately 50 µL) to the glass slide.
- Mix the erythrocytes with test serum well with a stick and spread to a circle (diameter 2 cm)
- By slightly rotating the slide, check for agglutination within 1 minute (reaction starts within seconds). Unspecific reactions might appear due to drying of the reaction-formation or if the slide is heated. Document the result.

### Tube Centrifugation Method

1. Use a 2% to 5% suspension of red blood cells in isotonic saline only (cells washed one time or up to three times with isotonic saline).
2. Add 100 µL (alternative: one drop = approximately 50 µL) of test serum to each tube.
3. Add 100 µL (alternative: one drop = approximately 50 µL) of appropriate cell suspension to each tube.
4. Mix well by slightly shaking.
5. Incubate tube at room temperature (+15 to +30 °C) for 15 min.
6. Centrifugation of tube for 1 minute at 2.000 rpm (approximately 800-1.000 g).
7. Gently resuspend the red cells and check macroscopically for agglutination within 3 minutes. Document the result.
8. Incubate all tubes with negative or weak result at +37 °C for 15 minutes.
9. Centrifugation of tube for 1 minute at 2.000 rpm (approximately 800-1.000 g).
10. Gently resuspend the red cells and check macroscopically for agglutination within 3 minutes. Document the result.
11. For all negative or weakly positive results, an indirect Coombs test must be performed.

### Indirect Coombs-test

12. Repeat step 1.-4. with fresh material and incubate tube at +37 °C for 15 min or continue after reading of result directly.
13. Wash red cells 3 times with (cold) isotonic saline.
14. Add 100 µL Anti-Human-Globulin-Serum (Coombs-serum) to each tube.
15. Centrifugation of tube for 1 min at 1.000 rpm (approximately 180-270 g).
16. Gently resuspend the red cells and check macroscopically for agglutination within 3 minutes. Document the result.

### INTERPRETATION OF RESULTS

"Slightly rotating / shaking" at Slide Method and at Tube Centrifugation Method

**Positive result (+):** visible agglutination of erythrocytes indicates the presence of one or more of the corresponding antigens.

**Negative result (-):** no visible agglutination of erythrocytes indicates the absence of all of the corresponding antigens.

### LIMITATIONS OF THE PROCEDURE

1. The "Procedures" and "Interpretation of Results" sections must be followed closely to assure the accuracy of the test result.
2. No valid conclusion concerning the test result can be reached, if controls with uncertain or false results occur.
3. Enzyme treated erythrocytes may react unspecific.
4. Due to variability of antigen expression, reactivity of this test serum against certain phenotypes may give weaker reactivity compared to control cells.
5. Red blood cells coated with alloantibodies or autoantibodies of the same or similar specificity as the test serum (i.e., cells that are positive in the direct antiglobulin-test (DAT)) may give weak reactions. In extreme cases, false-negative results may occur.
6. Cord blood is excluded from testing.
7. The reagent anti-CDE does not react with rG cells.
8. The IgM Anti-D component will react negatively with D category VI cells, the IgG Anti-D component reacts positive with these cells by the indirect antiglobulin method.

### PERFORMANCE CHARACTERISTICS

In compliance with Common Technical Specifications (CTS Commission Decision of 03. February 2009) a performance evaluation was conducted for the product. Different samples (donor, patient, panel blood) was used and compared with other products. The determined sensitivity and specificity was as follows

Total number of bleeding	Sensitivity	Total number of bleeding	Specificity
1.208	99,92%	705	99,86%

### LEGEND OF SYMBOLS

 Store from - to	 REF	Product Code	 Expiration Date	 Manufacturer as to 98/79/EU			
 LOT	Lot	 CLON	Clone(s)	 IVD	In vitro diagnostic medical device	 CE 0483	EU CE-symbol

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- Anti-C (Clone: MS-24 (IgM))
- Anti-D (Cloni: MS-201 (IgM), MS-26 (IgG))
- Anti-E (Clone: MS-80 (IgM))

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