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Anti-E (RH3), monoclonal IgM (human)

REF 690291A

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 Antitoxin GmbH, Industriestraße 88, 69245 Bammental

731-22-0515 Version 015 / 01. Sep 2014

In vitro diagnostic Medical Device
For Slide, Tube and Microplate Test and for combination

INSTRUCTIONS FOR USE

Anti-E (RH3), monoclonal IgM (human)

REF 690291A

INTENDED USE

Monoclonal agglutinating Anti-E reagent is produced from cell culture supernatants of heterohybridoma cell lines. The cells secrete an antibody of IgM-type, that reacts specifically with the corresponding antigen. The antibody is human protein. The reagent is used to determine whether red blood cells possess or lack the corresponding blood group antigens. The test sera are intended to be used by qualified technical personnel only.

PRINCIPLE OF PROCEDURE

The procedures used with this reagent are based on the principle of agglutination. Normal human erythrocytes, possessing the corresponding antigen, will agglutinate in the presence of the specific antibody directed towards the antigen.

REAGENT

The listed reagent is available with following clones:

Anti-E monoclonal, human IgM (clones: MS-258 and 906)

This reagent contains <0.1% (w/v) sodium azide as preservative. Additionally the reagent is comprised of active antibody, sodium chloride, macromolecules and bovine albumin.

CAUTION: Please handle reagent with proper care: This reagent was prepared from supernatants of cell cultures. Handle as if capable of transmitting infectious agents. The reagent contains sodium azide, which may be toxic and may react with lead or copper to form highly explosive salts. On disposal, flush with large quantities of water.

STORAGE REQUIREMENTS

Store at +2 to +8 °C (unopened / opened), or at room temperature (+15 to +30 °C) while in use. Do not use reagent beyond its labeled expiration date!

REMARKS

- Strength of positive reactions also depends on age of used blood.
- It is recommended that each lot of reagent be tested with appropriate positive and negative controls.
- Inappropriate storage impairs efficacy of the reagent.
- Overcentrifugation or undercentrifugation may lead to false results.
- Blood samples to be tested should be used as soon as possible. If a delay in testing occurs, samples should be stored at +2 to +8 °C. Blood drawn into sodium citrate or EDTA should be tested within 14 days.
- The procedures identified below are for manual testing only. When using automated or semi-automated instruments, follow the procedures that are contained in the operator's manual provided by the device manufacturer. Laboratories must follow approved validation procedures to demonstrate compatibility of this product on automated systems.
- For usage of this reagent all effective national laws, directives and guidelines must be followed.
[In Germany especially: „Richtlinien zur Gewinnung von Blut und Blutbestandteilen und zur Anwendung von Blutprodukten (Hämotherapie)“.]

REAGENT PREPARATION

There is no specific preparation of the reagent required. Use reagent directly from the vial.

PROCEDURE

Material required but not provided:

Slide Method

- Glass slide
- Pasteur pipette
- Mixing stick

Tube Centrifugation Method

- Test tubes, 10 x 75 mm or 12 x 75 mm
- Pipettes designed to deliver approximately 50 µL/100 µL
- Centrifuge
- Isotonic saline (0,85 - 0,9% Natriumchlorid)
- Disposable pipette tips

Microplate Method

- Microplate with U-bottom, (optional: pre-treated)
- Centrifuge, suitable for centrifugation of microplates
- Microplate-carrier for centrifuge (option)
- Microplate shaker (option)
- Microplatetest mirror (option)
- Pipette for 50 µL
- Isotonic saline (0.85 – 0.9 % sodium chloride)
- Disposable pipette tips

Test procedure

Slide Method

1. Use erythrocyte sediment or whole blood only.
2. Place one drop (approximately 50 µL) of reagent on a glass slide.
3. Using a Pasteur pipette add one drop erythrocyte sediment or whole blood (approximately 50 µL) to the glass slide.
4. Using a stick, mix the erythrocytes with reagent well and spread to a circle (diameter 2 cm).
5. By slightly rotating the slide, check for agglutination within 1 minute (reaction starts within seconds). Document results. Unspecific reactions might appear due to drying of the reaction-formation or if the slide is heated. Document the results.

Tube Centrifugation Method

1. Use a 2% to 5% suspension of red blood cells in isotonic saline only (cells washed one time or up to three times with isotonic saline).
2. Add 100 µL (alternative: one drop, approximately 50 µL) of reagent to each tube.
3. Add 100 µL (alternative: one drop, approximately 50 µL) of cell suspension to each tube.
4. Mix well by slightly shaking.
5. Incubate tube at room temperature (+15 to +30 °C) for 1-15 min.
6. Centrifugation of tube for 1 minute at 2.000 rpm (approximately 800 – 1.000 g).
7. Gently resuspend the red cells and within 3 minutes check macroscopically for agglutination. Document the result.

Microplate Method

1. Prepare a 2% to 5% suspension of red blood cells in isotonic saline. Cells can be washed one time or up to three times with isotonic saline.
2. Add one drop (50 µL) of reagent into each well.
3. Add one drop (50 µL) of cell suspension into each well.
4. Mix microplate on a microplate shaker.

Note: Recommended : (1) To mix: 10 to 30 seconds at medium strength. (2) To resuspend: 10 to 30 seconds at medium strength respectively with duration and strength, allowing a full resuspension of sediment without destruction of positive reactions.

5. Centrifugation of microplate. Recommended duration of centrifugation: 30 seconds at approximately 400 x g, or as long as necessary to produce best positive reactions with antigen positive cells, with used centrifuge and microplate, and still have an easy resuspension of antigen negative cells.

Note: The used centrifugal force should be the minimum force to produce a sediment and a clear supernatant. Because of multitude offered centrifuges a general speed or duration for best results can not be recommended. Each laboratory has to test and validate its equipment for speed and duration to obtain best results.

6. Loosen the erythrocyte knobs on the shaker and position microplate for reading.
7. Read microplate and document the result.
8. Tests with negative or doubtful results have to be incubated for 5 to 10 minutes at room temperature (+15 to +30 °C).
9. Repeat steps 5 to 7 after incubation at room temperature.

Note: All tests should be resuspended and evaluated directly after centrifugation.

INTERPRETATION OF RESULTS

"Slightly rotating / shaking " at Slide Method and at Tube Centrifugation Method, and using a microplate shaker at Microplate Method

- Positive result (+):** visible agglutination of erythrocytes indicates the presence of the corresponding antigen.
Negative result (-): no visible agglutination of erythrocytes indicates the absence of the corresponding antigen.

LIMITATIONS OF THE PROCEDURE

1. The "Procedure" and "Interpretation of Results" sections must be followed closely to assure the accuracy of the test results.
2. No valid conclusion concerning the test result can be reached, if controls with uncertain or false results occur.
3. Enzyme treated erythrocytes may react unspecifically.
4. Due to variability of antigen expression, reactivity of these reagents against certain phenotypes may produce a weaker reaction compared to control cells.
5. Red blood cells coated with alloantibodies or autoantibodies of the same or similar specificity as the reagent (i.e., cells that are positive in the direct antiglobulin-test [DAT]) may give weak reactions. In extreme cases, false-negative results may occur.

PERFORMANCE

In compliance with Common Technical Specifications (CTS Commission Decision of 03. February 2009) a performance evaluation was conducted. Different samples (donor, patient, panel blood) was used and compared with other products.

Product	Positive	False negative	Sensitivity	Negative	False positive	Specificity
Anti-E (clones: MS-258, 906)	327	0	100%	908	0	100%

LEGEND OF SYMBOLS

 Store from - to	 REF	Product Code	 Expiration Date	 Manufacturer as to 98/79/EU			
 LOT	Lot	 CLON	Clone(s)	 IVD	In vitro diagnostic medical device	 CE 0483	EU CE-symbol

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