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Antitoxin GmbH, Industriestraße 88, 69245 Bammental

730-22-0516 Version 016 / 30. Mrz 2015

## In vitro diagnostic Medical Device

For Slide, Tube and Microplate Test and for combination with Ortho BioVue® Reverse Diluent Cassettes

## INSTRUCTIONS FOR USE

**Anti-C (RH2), monoclonal IgM (human) OrthoClone**  
**Anti-E (RH3), monoclonal IgM (human) OrthoClone**  
**Anti-c (RH4), monoclonal IgM (human) OrthoClone**  
**Anti-e (RH5), monoclonal IgM (human) OrthoClone**

	690211A
	690321A
	690281A
	690361A

## INTENDED USE

Monoclonal agglutinating Anti-C, -E, -c, -e reagents are produced from cell culture supernatants of heterohybridoma cell lines. The cells secrete an antibody of IgM-type. The antibody in each case is human protein. The reagents are used to determine whether red blood cells possess or lack the corresponding blood group antigens C, E, c and e. The procedures described are for the listed manual methods and application with Ortho BioVue® Reverse Diluent cassettes on the ORTHO VISION™ Analyzer and ORTHO AutoVue® Innova. The test sera are intended to be used by qualified technical personnel only.

## PRINCIPLE OF PROCEDURE

The procedures used with these reagents are based on the principle of agglutination. Normal human erythrocytes, possessing the corresponding antigen, will agglutinate in the presence of the specific antibody directed towards the antigen.

## REAGENT

The listed reagents are available with following clones:

Anti-C monoclonal, human IgM (clone: MS-273)      Anti-c monoclonal, human IgM (clone: MS-35)  
Anti-E monoclonal, human IgM (clones: MS-12 und MS-260)      Anti-e monoclonal, human IgM (clones: MS-62 and MS-69)

These reagents contain <0.1% (w/v) sodium azide as preservative. Additionally these reagents are comprised of active antibody, sodium chloride, macromolecules and bovine albumin. CAUTION: Please handle reagent with proper care: These reagents were prepared from supernatants of cell cultures. Handle as if capable of transmitting infectious agents.

The reagents contain sodium azide, which may be toxic and may react with lead or copper to form highly explosive salts. On disposal, flush with large quantities of water.

## STORAGE REQUIREMENTS

Store at +2 to +8 °C (unopened / opened), or shortly at room temperature (+15 to +30 °C) while in use. If the devices ORTHO VISION™ Analyzer or ORTHO AutoVue® Innova are used, the reagents may remain a maximum of 5 days on the device (+15 bis +30 °C). Do not use reagents beyond the stated expiration date!

## REMARKS

- Strength of positive reactions also depends on age of used blood.
- It is recommended that each lot of reagents be tested with appropriate positive and negative controls.
- Inappropriate storage impairs efficacy of the reagent.
- Overcentrifugation or undercentrifugation may lead to false results.
- Blood samples to be tested should be used as soon as possible. If a delay in testing occurs, samples should be stored at +2 to +8 °C. Blood drawn into sodium citrate or EDTA should be tested within 14 days.
- The procedures described are for the listed manual methods and analyzers only. When using other automated or semi-automated instruments, follow the device manufacturer's procedures in the operator's manual. Laboratories must follow approved validation procedures.
- For usage of these test sera all effective national laws, directives and guidelines must be followed. [In Germany especially: „Richtlinien zur Gewinnung von Blut und Blutbestandteilen und zur Anwendung von Blutprodukten (Hämotherapie)“.]

## REAGENT PREPARATION

There is no specific preparation of the reagents required. Use reagents directly from the vial.

## PROCEDURE

Material required but not provided:

### Slide Method

- Glass slide
- Pasteur pipette
- Mixing stick

### Tube Centrifugation Method

- Test tubes, 10 x 75 mm or 12 x 75 mm
- Mikroliter pipette designed to deliver approximately 50µL/100 µL
- Centrifuge
- Isotonic saline (0,85 - 0,9% sodium chloride)
- Disposable pipette tips

### Microplate Method

- Microplate with U-bottom, (optional: pre-treated)
- Centrifuge, suitable for centrifugation of microplates
- Microplate-carrier for centrifuge (option)
- Microplate shaker
- Microplatetest mirror (option)
- Pipette for 50 µL
- Isotonic saline; (0,85 - 0,9 % sodium chloride)
- Disposable pipette tips

### Application of Ortho BioVue® Cassettes

- Ortho BioVue® System Centrifuge
- Isotonic saline (0,85 - 0,9 % sodium chloride)
- Mikroliter pipette for 10 µL, 40 µL and 50 µL
- Disposable pipette tips
- Ortho BioVue® System Work Rack
- Ortho BioVue® System Reverse Diluent Cassettes

## Test procedure

### Slide Method

- Use erythrocyte sediment or whole blood only.
- Place one drop (approximately 50 µL) of appropriate reagent on a glass slide.
- Using a Pasteur pipette add one drop erythrocyte sediment or whole blood (approximately 50 µL) to the glass slide.
- Using a stick, mix the erythrocytes with reagent well and spread to a circle (diameter 2 cm).
- By slightly rotating the slide, check for agglutination within 1 minute (reaction starts within seconds). Document results. Unspecific reactions might appear due to drying of the reaction-formation or if the slide is heated.

### Tube Centrifugation Method

- Use a 2% to 5% suspension of red blood cells in isotonic saline only (cells washed one time or up to three times with isotonic saline).
- Add 100 µL (alternative: one drop, approximately 50 µL) of appropriate reagent to each tube.
- Add 100 µL (alternative: one drop, approximately 50 µL) of appropriate cell suspension to each tube.

- Mix well by slightly shaking.
- Incubate tube at room temperature (+15 to +30 °C) for 1-15 min.
- Centrifugation of tube for 1 minute at 2.000 rpm (approximately 800-1.000 g).
- Gently resuspend the red cells and within 3 minutes check macroscopically for agglutination. Document the result.

#### Microplate Method

- Prepare a 2% to 5% suspension of red blood cells in isotonic saline. Cells can be washed up to three times with isotonic saline.
- Add one drop (approximately 50 µL) of appropriate reagent into each well.
- Add one drop (50 µL) of cell suspension into each well.
- Mix microplate on a microplate shaker.

**Note:** Recommended : (1) To mix: 10 to 30 seconds at medium strength. (2) To resuspend: 10 to 30 seconds at medium strength respectively with duration and strength, allowing a full resuspension of sediment without destruction of positive reactions.

- Centrifugation of microplate. Recommended duration of centrifugation: 30 seconds at approximately 400 x g, or as long as necessary to produce best positive reactions with antigen positive cells, with used centrifuge and microplate, and still have an easy resuspension of antigen negative cells.

**Note:** The used centrifugal force should be the minimum force to produce a sediment and a clear supernatant. Because of multitude offered centrifuges a general speed or duration for best results can not be recommended. Each laboratory has to test and validate its equipment for speed and duration to obtain best results.

- Loosen the erythrocyte knobs on the shaker and position microplate for reading.
- Read microplate and document the result.
- Tests with negative or doubtful results have to be incubated for 5 to 10 minutes at room temperature (+15 to +30 °C).
- Repeat steps 5 to 7 after incubation at room temperature.

**Note:** All tests should be resuspended and evaluated directly after centrifugation.

#### Manual Test Method in Ortho BioVue® Reverse Diluent Cassettes

The procedure described below is for manual Ortho BioVue® Cassettes testing only. Please refer to the instructions for use of the Ortho BioVue® Reverse Diluent Cassettes, sections STORAGE REQUIREMENT, PRECAUTIONS, SPECIMEN COLLECTION and PRÉPARATION, and INTERPRETATION OF RESULTS, as well as -where applicable- LIMITATIONS OF THE PROCEDURE.

- For preparation of the cassettes, please refer to the instructions for use for Ortho BioVue® Reverse Cassettes, chapter „procedure “ point 1 and 2.
- Prepare a 3 - 5% erythrocyte suspension in isotonic sodium chloride solution. Optionally, the erythrocytes may be washed up to three times in isotonic sodium chloride solution.
- Add 10 µL of a 3 to 5% suspension of test erythrocytes to the appropriate reaction chamber of the cassette.
- Add each 40 µL reagent to the appropriate reaction chambers of the cassette. The pipette tips shall not contact the rim of the reaction chambers. If this should adversely happen, please replace the tip with a new one, before you will fill the next reaction chamber.
- Centrifuge the cassette in the Ortho BioVue® System centrifuge.

**Note:** Centrifugation must be performed within 30 minutes after addition of specimens to the reaction chambers.

- After completion of test, read each column from the cassette's front and rear sides for agglutination and/or hemolysis.
- Protocol the final test result from the cassette side with the more clearly expressed positive reading.

#### Test in Ortho BioVue® Reverse Diluent Cassettes on ORTHO VISION™ Analyzer or ORTHO AutoVue® Innova

Please follow the instructions in the user manual of the respective analyzer.

#### INTERPRETATION OF RESULTS

"Slightly rotating / shaking" at Slide Method and at Tube Centrifugation Method, and using a microplate shaker at Microplate Method

- Positive result (+): visible agglutination of erythrocytes indicates the presence of the corresponding antigen.
- Negative result (-): no visible agglutination of erythrocytes indicates the absence of the corresponding antigens.

Cassettes: As for the manual method, for interpretation of results, please refer to the instructions for use of Ortho BioVue® Reverse Cassettes, chapter „INTERPRETATION OF RESULTS“. As for analyzers, reading and interpretation is effected automatically by each device.

#### LIMITATIONS OF THE PROCEDURE

- The "Procedure" and "Interpretation of Results" sections must be followed closely to assure the accuracy of the test results.
- No valid conclusion concerning the test result can be reached, if controls with uncertain or false results occur.
- Enzyme treated erythrocytes may react unspecifically.
- Due to variability of antigen expression, reactivity of these reagents against certain phenotypes may produce a weaker reaction compared to control cells.
- Red blood cells coated with alloantibodies or autoantibodies of the same or similar specificity as the reagent (i.e., cells that are positive in the direct antiglobulin-test [DAT]) may give weak reactions. In extreme cases, false-negative results may occur.
- Red blood cells coated with antibodies (cells that are positive in the direct antiglobulin-test (DAT)) may give false-positive results in cassettes method. These cells react positive without reagent too.
- Pay attention to all statements to limitations in inserts of used Ortho BioVue® Reverse Diluent Cassettes and the user manual of respective analyzer.
- Negative test results with anti-e (clones: MS-62, MS-69) in microplate method have to be confirmed by another technique, as false negative results may occur.
- With cord cells anti-c (clone: MS-33) may possess anti-I/i cold agglutinin activity and therefore show false positive results.


#### PERFORMANCE

In compliance with Common Technical Specifications (CTS Commission Decision of 03. February 2009) a performance evaluation was conducted. Different samples (donor, patient, panel blood) was used and compared with other products.

Product	Positive	False negative	Sensitivity	Negative	False positive	Specificity
Anti-C (clone: MS-273)	844	0	100%	391	0	100%
Anti-c (clone: MS-35)	1005	0	100%	230	0	100%
Anti-E (clones: MS-12, MS-260)	327	0	100%	908	0	100%
Anti-e (clones: MS-62, MS-69)	1417	0*	100%*	65	0	100%

\* for Anti-e (clones MS-62, MS-69) in microplate method 42 samples were tested false negative, so in this method sensitivity was 97,12% only. (see "Limitations of the procedure" item 8)

#### LEGEND OF SYMBOLS

	Store from - to		Product Code		Expiration Date		Manufacturer as to 98/79/EU
	Lot		Clone(s)		In vitro diagnostic medical device		EU CE-symbol

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**Anti-C (RH2), monoclonal IgM (human) OrthoClone**

**Anti-E (RH3), monoclonal IgM (human) OrthoClone**

**Anti-c (RH4), monoclonal IgM (human) OrthoClone**

**Anti-e (RH5), monoclonal IgM (human) OrthoClone**

	690211A
	690321A
	690281A
	690361A

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