

ABO-trol[®] forte plus

CE 0483

ImuMed
Lernfeldstraße 10
42699 Solingen



Verw. bis/
Exp. date:

Ch.-B.-Lot:

Lagern bei/
storage at:

ABO-trol[®] forte plus

30 gebrauchsfertige Bedside-Test-
karten zur ABO-Identitätsbestim-
mung vor der Bluttransfusion.

30 bedside-testcards ready for use
for verification of ABO-compatibi-
lity prior to blood transfusion.

ABO-trol / 730-01-0502



28440

30 gebrauchsfertige Bedside-Test-
karten zur ABO-Identitätsbestim-
mung vor der Bluttransfusion.
30 bedside-testcards ready for use
for verification of ABO-compatibi-
lity prior to blood transfusion.

30 bedside-testcards ready for use
for verification of ABO-compatibi-
lity prior to blood transfusion.
Pre-dipped fluid antisera with 70 µl
of each monoclonal anti-A (mouse),
anti-B (mouse) and anti-D (human)
in use for one determination.

Anti-A: Zelllinie celline BIRMA-1
Anti-D: Zelllinie celline RUM-1
Natriumazid sodium azide (<1 g/l)

ANTITOXIN GmbH
Industriestraße 88, 69245 Bammental (D)
FOR IN-VITRO-DIAGNOSTIC USE ONLY
+2 bis to +8 °C

522
30. 12. 2013